



Instant-View™  
rapid medical tests

## DRUG TEST - SALIVA



- Results in 5-7 Minutes
- Test Up to 9 Drugs at Once
- Accuracy Up to 99%
- Detects Parent THC at 50 ng/ml
- Minimal Sample Required
- Easy Collection Method
- Slim, Ergonomic Device Design

CE Marked  
Forensic Use Only in the USA  
US Patent D626249

Speed.  
Accuracy.  
Excellence.

Developed & Manufactured in the USA

US FDA Registered • ISO Certified



## Drug Tests

Up to 9 Analytes in One Test

### DRUG ANALYTE CHOICES

Product	Description	Sensitivity (ng/ml)
<b>1 STRIP SINGLE TEST</b>		
THC	Marijuana (parent)	50
<b>2 STRIPS IN MULTIPLE TESTS - 6-in-1</b>		
AMP	Amphetamine	50
MET	Methamphetamine	50
XTC	Ecstasy (MDMA)	50
COC	Cocaine	20
OPI	Opiates	40
THC	Marijuana (parent)	50
<b>2 STRIPS IN MULTIPLE TESTS - 7-in-1</b>		
AMP	Amphetamine	50
MET	Methamphetamine	50
XTC	Ecstasy (MDMA)	50
BZD	Benzodiazepine	20
COC	Cocaine	20
OPI	Opiates	40
THC	Marijuana (parent)	50
<b>2 STRIPS IN MULTIPLE TESTS - 9-in-1</b>		
AMP	Amphetamine	50
MET	Methamphetamine	50
XTC	Ecstasy (MDMA)	50
BZD	Benzodiazepine	20
COC	Cocaine	20
OPI	Opiates	40
MTD	Methadone	20
PCP	Phencyclidine	10
THC	Marijuana (parent)	50

### Product Format

- Saliva Card

### Product Registration Status

- **CE** Marked
- Forensic / Research

Made in  
the USA